

SPECIAL PRODUCTS FOR OEDEMA THERAPY

Compression bandaging with Juzo[®] Softcompress



Juzo[®]. Freedom in Motion

Our customer service department – always there for you

Service with a capital S

To us, customer service means being reachable and always having a sympathetic ear. Our service teams will deal with your concerns with passion, flexibility and reliability. You can reach us by phone **MON to FRI** or 24/7 by email or fax.

Customer service also means expert advice. Take advantage of our more than 100 years of experience in compression therapy. Our friendly and dedicated staff will be happy to assist you with all your enquiries.

We look forward to your call!

How to reach us

- General information
- Product advice
- Order acceptance
- Order status
- Alterations and repairs

Customer service



Edith Männer
Customer Service
Manager International

How to reach us

Juzo Customer Service International	
Office hours	7:00 am - 6:00 pm
Phone	+49 82 51 90 13 00
Fax	+49 82 51 90 13 21
E-mail	export@juzo.de

Juzo®

More than 100 years of experience
in compression therapy

Summary



Compression bandaging

Material and characteristics

04

How it works

05



Bandaging aids

Uni sizes

06

Standard sizes

06

Custom-made items

06

Compression dressings

09

Compression sheets

09



Lymph and bandaging sets

Lymph set

10

Bandaging set

11



Accessories

Breast pad and bandage

12

Genital pad

12

Joint and malleo pad

13

Hand pad

13



Material and characteristics

Using Juzo Softcompress can significantly improve the treatment outcome. Compared to the often used foam or synthetic cotton, Juzo Softcompress is characterized by its particular skin friendliness and durability. All Juzo Softcompress Bandaging Aids and Accessories are manufactured in Germany.

Characteristics

- Especially skin-friendly as the surface consists of 100 % cotton tricot
- Core made from foam
- Vertical seams
- Washable up to 60 °C
- Breathable
- Saves time when applying bandaging
- Makes self-bandaging easier
- Provides long-term cost savings

The following wash instructions apply to Juzo Softcompress compression material:



How it works

These comfortable compression aids can be worn under short-stretch bandages or compression stockings, and provide additional compression, cushioning and an ideal pressure distribution.

Juzo Softcompress Bandaging Aids are employed as an effective treatment aid for lymphedema and lipedema. The vertically-stitched foam channels create differences in pressure, thereby promoting lymph drainage.

The bandaging aids avoid any constriction by the edges of the compression bandage and protect the bends of the joints, while at the same time preventing soft tissue fibrosis. These characteristics and their effect are what make Juzo Softcompress so unique.



Visible effect of the Juzo Softcompress bandaging aid on a secondary leg lymphedema





Bandaging aids

The comfortable Juzo Softcompress Bandaging Aids can be worn under short-stretch bandages or compression stockings, and provide additional compression, cushioning and an ideal pressure distribution.



They are available in uni sizes, standard sizes and custom-made items for fast and effective bandaging. They make self-bandaging easy for the patient.

Custom-made bandages are ideal for night time therapy.

Juzo Softcompress Bandaging Aid Leg

	Art. no.	Description	Uni size (length cannot be cut)	Standard size 1-5 (length can be cut)	Custom-made
	6830	Bandaging aid lower leg incl. 2 malleo pads		●	●
	6831	Bandaging aid lower leg incl. 2 malleo pads	●		
	6820	Bandaging aid lower leg and thigh			●
	6821	Bandaging aid lower leg and thigh incl. 2 malleo pads, 1 joint pad, 1 compression roll (10 cm x 4 m) and 1 compression roll (12 cm x 4 m)	●		

Juzo Softcompress Bandaging Aid Arm

	Art. no.	Description	Uni size (length can be cut)	Standard size 1-6 (length can be cut)	Custom-made
	6822	Bandaging aid arm incl. 1 hand pad and 1 joint pad		●	●
	6823	Bandaging aid arm incl. 1 hand pad and 1 joint pad	●		

Size charts for standard sizes see page 9.
For custom-made items please use the Juzo Order Forms.



Bandaging aid lower leg
in standard size



Bandaging aid lower leg and thigh
in uni size



Bandaging aid lower leg and thigh
in custom-made



Bandaging aid arm
in custom-made




Bandaging aid arm
in uni size

Size charts

Bandaging aid lower leg in standard sizes


Größe	cA	cY	cB	cB1	cC	cD
1	18,5-20	28,30,5	18,5-20	23,5-25,5	29-31,5	27,29,5
2	20-22	29-31,5	20-22	25,5-28	31,5-34,5	29,5-32,5
3	22-24	31-34	22-24	27,5-30	34,5-37,5	32,5-35,5
4	24-26	33-36	24-26	30-33	37,5-41	35,5-39
5	26-28,5	35-38,5	26-28,5	32,5-35,5	40,5-44,5	38,5-42



The diagram shows a side view of a lower leg. Horizontal lines indicate the measurement points for the bandaging aid: cA at the ankle, cY at the heel, cB at the medial malleolus, cB1 at the lateral malleolus, cC at the medial malleolus, and cD at the lateral malleolus.

Bandaging aid arm in standard sizes

	Umfang					Länge	
Größe	cG	cF	cE	cD	cC	lcG	lcE
1	30-33	28-30,5	25-27,5	23-25	14,5-16	46-51	21-25
2	32-35	30-33	27-29,5	24,5-27	16-17,5	46-51	21-25
3	34-37	32-35	29-31,5	26-28,5	17,5-19	46-51	21-25
4	36-39,5	34-37	30,5-33,5	27,5-30	19-20,5	46-51	21-25
5	38-41,5	36-39,5	32-35	29-31,5	20,5-22	46-51	21-25
6	40-44	38-41,5	33,5-36,5	30,5-33,5	22-24	46-51	21-25



The diagram shows a side view of an arm. Horizontal lines indicate the measurement points for the bandaging aid: cG at the shoulder, cF at the elbow, cE at the wrist, cD at the wrist, cC at the wrist, lcG at the shoulder, and lcE at the wrist.

Juzo Softcompress Compression Roll

These thin, very soft compression dressings are ideal for use under lymphological and phlebological compression bandages. The Softcompress Compression Rolls distribute the pressure of the short-stretch bandages, thereby avoiding constriction.

Size	Thickness	Dimension	Art. no.
1	5 mm	8 cm x 4 m	6835
2	5 mm	10 cm x 4 m	6835
3	5 mm	12 cm x 4 m	6835
4	5 mm	16 cm x 4 m	6835



Juzo Softcompress Compression Sheet

Juzo Softcompress Compression Sheets can be cut to size for use in various body regions.

Size	Dimension	Art. no.
1	65 cm x 65 cm	6836
2	125 cm x 100 cm	6836



Lymph and bandaging set

The Juzo Softcompress Bandaging Sets and Lymph Sets offer a complete solution for lymphological compression bandaging. Unlike conventional sets, they also include the innovative Juzo Softcompress Compression Aids. The sets are suitable for compression bandaging at home. Using these compression aids can significantly improve the treatment outcome.

The lymph sets come with the ready-to-use bandaging aids, while the more economical bandage sets contain the Softcompress Compression Rolls.

Juzo Softcompress Lymph Set

Lymph set leg small

Art. no. 6800 KL

1 bandaging aid lower leg in uni size,
20 gauze bandages 4 m (for compression of toes),
2 short-stretch bandages 6 cm x 5 m,
2 short-stretch bandages 8 cm x 5 m,
2 short-stretch bandages 10 cm x 5 m,
4 short-stretch bandages 12 cm x 5 m and
2 malleo pads



Lymph set leg large

Art. no. 6800 GR

1 bandaging aid lower leg in uni size,
1 compression roll 10 cm x 4 m,
1 Kompressionsbinde 12 cm x 4 m,
20 gauze bandages 4 m (for compression of toes),
2 short-stretch bandages 6 cm x 5 m,
2 short-stretch bandages 8 cm x 5 m,
4 short-stretch bandages 10 cm x 5 m,
6 short-stretch bandages 12 cm x 5 m,
2 malleo pads and
1 joint pad



Lymph set arm

Art. no. 6801

1 bandaging aid arm in uni size,
20 gauze bandages 6 cm (for compression of fingers),
2 short-stretch bandages 6 cm x 5 m
2 short-stretch bandages 10 cm x 5 m
3 short-stretch bandages 12 cm x 5 m
1 hand pad and
1 joint pad



Juzo Softcompress Bandaging Set

Bandaging set leg small

Art. no. 6810 KL

1 compression roll 10 cm x 4 m,
20 gauze bandages 4 m (for compression of toes),
1 short-stretch bandages 6 cm x 5 m,
1 short-stretch bandages 8 cm x 5 m,
3 short-stretch bandages 10 cm x 5 m and
2 malleo pads



Bandaging set leg large

Art. no. 6810 GR

1 compression roll 10 cm x 4 m,
1 compression roll 12 cm x 4 m,
20 gauze bandages 4 m (for compression of toes),
1 short-stretch bandages 6 cm x 5 m,
1 short-stretch bandages 8 cm x 5 m,
3 short-stretch bandages 10 cm x 5 m,
4 short-stretch bandages 12 cm x 5 m,
2 malleo pads and
1 joint pad



Bandaging set arm

Art. no. 6811

1 compression roll 10 cm x 4 m,
20 gauze bandages 6 cm (for compression of fingers),
1 short-stretch bandages 6 cm x 5 m
1 short-stretch bandages 10 cm x 5 m
1 short-stretch bandages 12 cm x 5 m
1 hand pad and
1 joint pad



The gauze bandages and short-stretch bandages in the lymph sets and bandage sets are quality products from the company Karl-Otto-Braun (KOB). They are also available separately.


Accessories

The Juzo Softcompress Accessories are special pads that can be worn under compression bandages. They are especially suitable for treating concave areas and particularly localized edemas.

Juzo Softcompress Breast Pad

The Juzo Softcompress Breast Pad serves as a comfortable liner for a compression bra or chest bandage, e.g. for breast or chest wall edemas. It is available in three standard sizes and in custom-made.

Size	Cup size	Art. no.
1	A-B	6837
2	C-D	6837
3	E-F	6837



Juzo Softcompress Breast Bandage

The Juzo Softcompress breast bandage combines the positive characteristics of the Juzo Softcompress material with those of a compression bra. The breast bandage can, for example, be used for breast or chest wall edemas. It is available in five standard sizes and in custom-made.

Size	Cup size	Art. no.
1	A	6838
2	B	6838
3	C	6838
4	D	6838
5	E	6838
6	F	6838



Juzo Softcompress Genital Pad for Women

The Juzo Softcompress Pad for the female genital area is worn underneath the compression stocking. It provides a targeted increase in pressure, cushioning and an ideal pressure distribution when treating genital edemas. It is available in four standard sizes.

Size	Dimensions	Art. no.
1	20 cm x 20 cm	6839
2	30 cm x 28 cm	6839
3	40 cm x 32 cm	6839



Juzo Softcompress Genital Pad for Men

Art. no. 6840

The Juzo Softcompress Pad for the male genital area is worn underneath the compression stocking. It provides a targeted increase in pressure, cushioning and an ideal pressure distribution when treating genital edemas. The Juzo Softcompress genital pad for men is only available in custom-made.



Juzo Softcompress Joint Pad

Art. no. 6842

Juzo Softcompress Joint Pads cushion and soften the area in the bend of the joint, for example at the knee or elbow. These liners can also be used for various other applications.



Juzo Softcompress Malleo Pad

Art. no. 6843

Juzo Softcompress Malleo Pads cushion the Bisgaard (perimalleolar) region of the ankle. The pads are delivered pairwise.



Juzo Softcompress Hand Pad

Art. no. 6844

Juzo Softcompress Hand Pads are especially for the back of the hand with a two finger openings.



Julius Zorn GmbH

Juliusplatz 1
86551 Aichach
Germany
info@juzo.com

www.juzo.com

