

Pituusmitat, cm

Ympärysmitat, cm



		Vakio koot						
		I XS	II S	III M	IV L	V XL	VI XXL	
<b>Kompressio hihat</b> Art. 2001, 2002 3511, 3512		cG max	30-33	32-35	34-37	36-39,5	38-41,5	40-44
		cG	24-26	26-28,5	28-30,5	30-33	32-35	34-37
		cF max	28-30,5	30-33	32-35	34-37	36-39,5	38-41,5
		cF	23-25	24,5-27	26-28,5	27,5-30	29-31,5	30,5-33,5
		cE max	25-27,5	27-29,5	29-31,5	30,5-33,5	32-35	33,5-36,5
		cE	23-25	24,5-27	26-28,5	27,5-30	29-31,5	30,5-33,5
cD max	25-27,5	27-29,5	29-31,5	30,5-33,5	32-35	33,5-36,5		
cD	23-25	24,5-27	26-28,5	27,5-30	29-31,5	30,5-33,5		
cC' cC	14,5-16	16-17,5	17,5-19	19-20,5	20,5-22	22-24		
<b>JuzoFlex® Kynnärpäätuki</b> Art. 3022		cF	24	26	28	30	32	34
		cE	23	24,5	26	27,5	29	30,5
		cD	23	24,5	26	27,5	29	30,5
<b>JuzoFlex® Rannetuki</b> Art. 3022		cC	14,5	16	17,5	19	20,5	22

Julius Zorn GmbH  
 Postfach 1280, 86543 Aichach/Saksa  
 Juliusplatz 1, 86551 Aichach/Saksa

Puh. 00 49/82 51/90 13 00  
 Fax 00 49/82 51/90 13 21

E-mail: info@juzo.de  
 www.juzo.com

Ab Napra-Rehab Oy  
 Kuskgränden 3  
 64200 Närpiö/Finland

Puh 06-2 24 36 22  
 Fax 06-2 24 36 62

E-mail: info@napra-rehab.fi  
 www.napra-rehab.fi