

magnide on/off

simply smile

NEW

Magnide® On/Off

The ultimate application aid for putting on and taking off closed toe compression stockings and tights in one product.



Problem definition

For both caregivers and patients the problem of putting on and taking off compression stockings and tights is one which occurs on a daily basis. For caregivers this is also a major source of physical strain. Even patients who are able to carry out these two proceedings independently may find the process an increasingly arduous one. This often has a direct effect on therapy compliance.

Solution

The Magnide On/Off has been specifically designed to alleviate the problems experienced when putting on and taking off compression stockings and tights with closed toe. This product is made of smooth material, which has been treated with a special coating. These qualities create a smooth aid that significantly reduces friction when putting on and taking off stockings.

The friction reduction of the Magnide On/Off helps to minimise the physical strain otherwise experienced by caregivers and patients. In turn, this has a positive effect on the therapy compliance of patients. With clear and simple instructions-for-use patients may long enjoy the freedom and convenience of independently putting on and taking off compression hosiery.

Users

The Magnide ON/Off provides an effective solution for patients who wear compression stockings and tights with a closed toe. Thanks to the Magnide On/Off hosiery can now be put on and taken off without difficulties. Furthermore, the product is an essential aid for caregivers who are faced with the prospect of helping their patients putting on and taking off compression hosiery on a daily basis.


Magnide® On/Off

The ultimate solution in one product

Product benefits

- Two in one
- Top quality - outstanding design
- Reduction of physical strain for both caregiver and patient
- Increase in patient's independence and therapy compliance
- Stockings retain therapeutic graduated compression longer
- Ease of use
- Saves time and money
- Durable
- Easy to carry

Product range





Art. nr.	size	shoe size	
E01530	Medium	≤ 40	XL also suitable for extra wide feet with cA upto 28 cm 
E01531	Large	≥ 41 t/m 45	
E01532	Extra Large	≥ 46	




Transversely installed magnetic pocket



Nested eyelet for elastic cord with puller



- 1 Place the Magnide On/Off across the foot. Ensure that the orange flap near the toes is well folded back, so the magnet makes contact at the bottom.
- 2 Pull the stocking over the Magnide On/Off. Advice: don't roll up the stocking because that makes it even more difficult to put on!
- 3 Slide the stocking over the Magnide On/Off until the heel of the stocking is in the right position.
- 4 Make sure that the light green loop sticks out above the stocking so it can be grabbed.
- 5 Wiggle your toes. Pull the light green loop upwards to remove the Magnide On/Off from the stocking. Regrab if possible when removing the Magnide On/Off.



- 1 Slide the Magnide On/Off over the stocking.
- 2 Firmly pull the orange tab so that the Magnide On/Off is allocated just above the ankle.
- 3 Fold the stocking over the Magnide On/Off as far as possible.
- 4 Starting at the heel, gradually slide the stocking and the Magnide On/Off from the foot.