



**Juzo**

Freedom in Motion

# Lymphology Maintenance therapy whilst sleeping

Securing treatment success at night with Juzo SleepCompress

Lymphology

# Juzo SleepCompress

## Avoid disrupting your treatment

During the maintenance phase of oedema therapy, Juzo SleepCompress is the perfect addition to a flat knit compression garment. As soon as the compression garment has been taken off, it is important to make sure that the oedema does not increase again. At night and during rest phases, Juzo SleepCompress helps to secure the treatment's success, thus making oedema therapy even more effective and comfortable for patients.



### Foam channels

follow lymphatic pathways and promote lymphatic drainage

### Anatomical angle

reduces creasing and prevents pressure points as well as skin irritation

### Removable touch fasteners

ensure easy and quick donning

### Cooling material

cools the skin and reduces sweating at night



Cared for  
around the clock



For rest  
phases



Comfortable  
and  
convenient



**Thumb stub**  
prevents the oedema  
from shifting



## Can be used anywhere, at any time

Juzo SleepCompress is part of user-friendly self-care in addition to oedema therapy. The product is perfect for use at night and during other rest phases when garments are often omitted for reasons of convenience. The Juzo SleepCompress is easy and quick to put on; it is convenient as well as comfortable, and can be taken anywhere. In this way, it makes self-management and compliance easier for patients – both at home and on the go, as well as when travelling.

### Successful self-management

- lining and touch fasteners make it easier to put on
- individually adjustable; can be adapted at any time
- adapts to fluctuations in circumference that occur over the course of treatment




### Lining

helps with securing when the garment is put on

## The perfect good-night garment

The night is for recharging your batteries so you can start the day fresh and strong – Juzo SleepCompress makes it happen. The product is easy to use and exceptionally comfortable to wear. Juzo SleepCompress also reduces feelings of tension and pressure caused by oedema. This means you will feel protected and perfectly cared for even during rest phases.


## Use of Juzo lymphology products during complex manual decongestive therapy





## Everything going the right way

The structure of Juzo SleepCompress models itself on the lymphatic pathways. The stitched foam channels guide and support lymphatic drainage, thus ensuring decongestion. This reduces swelling, meaning the pressure caused by oedemas and the feeling of tension are decreased.



### Supporting the lymphatic system

The foam channels ensure the perfect pressure distribution and support lymphatic drainage.



#### Indications, contraindications and material composition

See usage and care instructions at [juzo.com/sleepcompress](https://www.juzo.com/sleepcompress)

#### Colour



Black



Washable  
60 degrees



Suitable for  
tumble-drying



# Juzo SleepCompress · Lower leg, leg

## Size chart

### Juzo SleepCompress · Lower leg, leg

Art. 6873 · Art. 6874		Circumference c Length P	Standard sizes · Measurements in cm				
			1/XS	2/S	3/M	4/L	5/XL
<b>PAG AG</b> long 72–83 normal 66–72		<b>cG</b>	40–46	46–52	52–58	58–64	64–70
		<b>cF</b>	38–44	44–50	50–56	56–62	62–68
		<b>cE</b>	28–33	33–38	38–43	43–48	48–53
<b>PAD AD</b> long 40–46 normal 34–40		<b>cD</b>	25–29	29–33	33–37	37–42	42–47
		<b>cC</b>	27–31	31–35	35–39	39–44	44–49
		<b>cB1</b>	21.5–25	25–28.5	28.5–32	32–36	36–40
		<b>cB</b>	17–20	20–23	23–26	26–29.5	29.5–32
		<b>cY</b>	26–29	29–32	32–35	35–38.5	38.5–42
		<b>cA</b>	17–20	20–23	23–26	26–29.5	29.5–32

## Product Range

	Art.	Model	Standard sizes
	<b>6873</b>	<b>Juzo SleepCompress Lower leg</b>	Each
	normal 34–40 long 40–46	left right	Size 1
			Size 2
			Size 3
			Size 4
		Size 5	
	<b>6874</b>	<b>Juzo SleepCompress Leg</b>	Each
	normal 66–72 long 72–83	left right	Size 1
			Size 2
			Size 3
			Size 4
		Size 5	
	<b>9572</b>	<b>Touch fastener pieces</b>	Set: 6 Each
	Size: 8.5 x 7.5 cm	Lower leg	Uni size
	<b>9573</b>	<b>Touch fastener pieces</b>	Set: 9 Each
	Size: 8.5 x 7.5 cm	Leg	Uni size

# Juzo SleepCompress · Arm

## Size chart

### Juzo SleepCompress · Arm

Art. 6872	Circumference c Length ℓ	Standard sizes · Measurements in cm				
		1/XS	2/S	3/M	4/L	5/XL
<p>PCG long 47–51 normal 42–47 short 36–42</p>	cG	28–32	31–35	34–38	37–41	40–44
	cF	26–29	29–32	32–35	35–38	38–41
	cE	23–26	26–29	29–32	32–35	35–38
	cD	20–23	23–26	26–29	29–32	32–35
	cC	12.5–15	15–17.5	17.5–20	20–22.5	22.5–25

## Product Range

	Art.	Model	Standard sizes
	<b>6872</b>	<b>Juzo SleepCompress Arm</b>	Each
	short 36–42 normal 42–47 long 47–51	left right	Size 1
			Size 2
			Size 3
			Size 4
	<b>9572</b>	<b>Touch fastener pieces</b>	Set: 6 Each
	Size: 8.5 x 7.5 cm	Arm	Uni size



Arm · Art. 6872




Lower leg · Art. 6873



Leg · Art. 6874







Further information is available at  
[juzo.com/sleepcompress](https://juzo.com/sleepcompress)

### Juzo Customer Service International

Phone: +49 (0)8251 901 300 · E-mail: [order@juzo.com](mailto:order@juzo.com) · Shop: [portal.juzo.com](https://portal.juzo.com)



Julius Zorn GmbH · Juliusplatz 1 · 86551 Aichach · Germany · [juzo.com](https://juzo.com)

ENG - 14999003/3/2 LY - 11/2024  
Subject to change, no liability for errors.

